



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
WATER PUMP PLIER - (GROOVE JOINT TYPE)

**Product Code :**  
LBNY-0031-2200010



**Description :**

Induction Hardened Jaws.  
Ribbed Handle for Slip Free Gripping.  
Drop Forged from Selected Carbon Steel.  
Available in With or Without Dip Sleeve.

**Technical Specification :**

SIZE (INCH): 8, 10, 12  
SIZE (MM.): 200, 250, 300



**Laboratory Instrument India**