



Email : sales@laboratoryinstrumentindia.com

Product Name :
JOINT WATER PUMP PLIER (CHANNEL TYPE)

Product Code :
LBNY-0031-220011



Description :

Drop Forged from Selected Carbon Steel.
Ribbed Handle for Slip Free Gripping.
Induction Hardened Jaws.
Available in Powder Coated and Dip Insulation Sleeve also.

Technical Specification :

SIZE (INCH): 8, 10, 12
SIZE (MM.): 200, 250, 300



Laboratory Instrument India