



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
LOCK GRIP PLIER

**Product Code :**  
LBNY-0031-2200012



**Description :**

Induction Hardened and Tempered Jaws.  
Jaws made from Selected Carbon Steel.  
Available in Chrome Plated / Powder Black Coated.

**Technical Specification :**

SIZE (INCH): 5, 7, 10  
SIZE (MM.): 125, 175, 250



**Laboratory Instrument India**