



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Automatic Precipitation (Tipping Bucket) to be Installed At  
The River Gauging Stations

**Product Code :**  
LBNY-0002-830009



**Description :**

Automatic Precipitation (Tipping Bucket) to be Installed At The River Gauging Stations

**Technical Specification :**

Automatic Precipitation (tipping Bucket). To be installed at the river gauging stations:  
high-quality tipping bucket rain gauge for measuring rainfall and precipitation in urban and rural locations. Due to the integrated syphon, the gauge delivers high levels of accuracy across a broad range of rainfall intensities. The Tipping Bucket tried and proven design ensures long-term, accurate and repeatable results. It is manufactured from high quality, durable materials ensuring long-term stability in the harshest of environments. It consists of a robust powder-coated aluminium enclosure, an aluminium base, and stainless steel ?nger ?lter and fasteners.



**Laboratory Instrument India**