



Email : sales@laboratoryinstrumentindia.com

Product Name :
Radar sensor Telemetry System

Product Code :
LBNY-0002-830010



Description :

Radar sensor Telemetry System

Technical Specification :

Radar sensor (Telemetry System)

Cable (Total length in m)

50m -No. 1

90m – No. 2

High performance

Measurements are unaffected by air temperature, humidity, flood events, floating debris, or contaminated water; reduces the likelihood of missing data and reduces data post processing

Low profile design

Unobtrusive appearance is ideal for urban installation sites or sites prone to vandalism.

Low maintenance

Flat antenna design eliminates nesting areas for insects and periodic maintenance requirement.

Simple installation

Light weight compact design facilitates easy mounting on bridges, extension mounting bracket or inside a small protective housing.

Low power consumption

Ideal for remote or solar powered sites; requires only 12 mA when active.



Laboratory Instrument India