



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
OIL CANS PERFETTO TYPE

**Product Code :**  
LBNY-0035-220001



**Description :**

A dome shaped body from heavy gauge has a modern design Detachable rigid & flex spouts Body is seamless drawn.

**Technical Specification :**

CAPACITY CC	CAPACITY OZ
200	6
300	10
400	13
500	16
700	23
1000	1 quart



**Laboratory Instrument India**