



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
OIL CANS PINT TYPE

**Product Code :**  
LBNY-0035-220002



**Description :**

Body of oil cans drawn from coated sheet heavy gauge.  
Detachable rigid & flex spouts.  
Body is seamless drawn.

**Technical Specification :**

**CAPACITY**

1/4 pint  
1/2 pint  
3/4 pint  
1 pint  
1 ½ pint  
2 pint



**Laboratory Instrument India**