



Email : sales@laboratoryinstrumentindia.com

Product Name :
Hydrostatic Pressure sensor Telemetry System

Product Code :
LBNY-0002-830011



Description :

Hydrostatic Pressure sensor Telemetry System

Technical Specification :

Hydrostatic Pressure sensor (Telemetry System)

Hydrostatic Pressure Sensor is a programmable submersible level transmitter for measuring water level and temperature in groundwater and surface water. It is low power, robust, and especially designed for use in harsh environments, ensuring longterm measurement stability.

The Hydrostatic Pressure Sensor provides very accurate results: The sensor is vented and automatically self-compensates for barometric pressure changes, i.e. it eliminates errors up to several inches that can be caused by pressure changes. Additionally temperature compensation is calculated from the measured water temperature. The fast reacting membrane technology of the Hydrostatic Pressure Sensor provides long-term measurement stability.

Cable (different lengths in m)

70m -No. 1

60m – No. 3



Laboratory Instrument India

