



Email : sales@laboratoryinstrumentindia.com

Product Name :
HAND OPERATED PUMP

Product Code :
LBNY-0035-220006



Description :

Lever Operated pump mounted on cover of drum.
Comes complete with 2-2 meter 0 high pressure hose fitted with accessories.
Working pressure up to 4000 PSI.
Pump has steel constructed cast head and solid steel pump chamber.
Quick & easy openable cover
Removable handle.

Technical Specification :

CAPACITY KGS:
5
10
15
20.



Laboratory Instrument India

