



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
ROTARY BARREL PUMP

**Product Code :**  
LBNY-0035-220007



**Description :**

The Rotary Pumps is superior to any other with features such as Dual Directional Operation which allows the pump to both empty, as well as refill containers Pump body is study cast iron with precisely machined cast iron vanes for better draw & smoother delivery Pump seals are made from Hi Nitrile rubber fits 15 gallons (50 litres) to 55 gallons (205 liter) bung nut & steel discharge spout. Most popularly used pump worldwide.

**Technical Specification :**

SIZE	DESCRIPTION	SUCTION TUBE
3/4"	STD. Rigid	
1"	STD. Rigid	
1"	STD. with flexible hose 1 mtr.	Telescopic pipe
1"	STD without flexible hose	Telescopic pipe
1 1/2"	STD Rigid	



**Laboratory Instrument India**