



Email : sales@laboratoryinstrumentindia.com

Product Name :
Main Cable For Cableway System Galvanized Tow cable
10 mm

Product Code :
LBNY-0002-830012



Description :

Main Cable For Cableway System Galvanized Tow cable 10 mm

Technical Specification :

Main cable for cableway system (Galvanized Tow cable 10 mm)

Material: Premium steel wire car tow rope, performs stronger and provides a more consistent pulling force than any other material. The looped ends are well built and highly strengthened to take the tension at the connection points.

The blue PVC hose attaches to the surface of the steel wire, effectively protect from scratches, beautiful and durable.

Firm and Solid Lock Hook: Large forging hook with lock, more firmer and safer, To prevent the hook off the trailer.

This item is a has tool for drivers when they get stuck in the mud, sand or snow and other emergency.



Laboratory Instrument India

