



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Traversing Cable for Cableway System Galvanized Tow  
Cable 8mm

**Product Code :**  
LBNY-0002-830013



**Description :**

Traversing Cable for Cableway System Galvanized Tow Cable 8mm

**Technical Specification :**

Traversing cable for cableway system (Galvanized Tow cable 8mm)  
Steel wire rope with 2 Drop Forged Steel Hooks, Strong & durable Galvanization technology makes the rope and hooks not easy to rust, increase its service life, It is a much necessary tool for drivers when cars get stuck in the mud sand or snow in Emergency, Heavy duty steel rope as long as 4 meters can withstand weight of 5 tons car, firm enough for use.



**Laboratory Instrument India**