



Email : sales@laboratoryinstrumentindia.com

Product Name :
DROP FORGEL BENCH VICE FIXED/SWIVEL

Product Code :
LBNY-0037-220001



Description :

All Parts are Forged From Selected Material
Corrosion Resistant Nickel Plated Handle
Hardened Fixed Serrated steel Jaws / Replacable jaws
Hardness of 42-50 HRC
Pipe vice jaw, large Hardened Anvil
Lead Screw with Double thread so as to increase the Efficiency
360 degree Swivel Base also available.
Parts are machined on VMC for Accurate & Precision Working
Nut in Housed in body and gives long service & can be replaced.

Technical Specification :

SIZE INCH	MAX. JAW OPENING	WEIGHT KGS. (FIXED)	WEIGHT KGS. (SWIVEL)
100	100	6.5	7.5
125	125	11.5	13.75
150	150	21.0	22.5



Laboratory Instrument India

