




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
BABY VICE CLAMP TYPE / SWIVEL BASE

**Product Code :**  
LBNY-0037-220003



**Description :**

Made out of Cast Iron with Guide Rods

**Technical Specification :**

SIZE FIXED	SIZE SWIVEL	SIZE MM
1	25	
1.1/2	1.1/2	40
2	2	50
2.1/2	2.1/2	60
3	3	70



**Laboratory Instrument India**