



Email : sales@laboratoryinstrumentindia.com

Product Name :
CHAIN PIPE VICE

Product Code :
LBNY-0037-220007



Description :

Compact, portable and easily fixed to bench, post stand.
Cast Iron Body.
High tensile steel chain with positive adjustment nut and handle.
Replacement drop forged steel jaws.
Chain vise that could be bolted to the back of a pickup but also taken off and bolted to a tripod stand.
Used for powerful Grip.

Technical Specification :

SIZE INCH	SIZE MM
2	50
3	75
4	100
6	150
8	200



Laboratory Instrument India

