



Email : sales@laboratoryinstrumentindia.com

Product Name :
Data Collection Platform

Product Code :
LBNY-0002-830015



Description :

Data Collection Platform

Technical Specification :

Data Collection Platform

Datalogger, memory card, internal battery, and spare parts

Enclosure and accessories .

Built-in 3.3V regulator is a reliable reference voltage, and a reliable operation, requires a lot of power to run the SD card.

Welded joints, circuits or sensors prototyping area.

Real-time clock (RTC) to maintain the time, even when the arduino is unplugged.

SD card interface with formatted card. 3.3V level shifter circuit to prevent damage to your SD card.



Laboratory Instrument India