




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Centripetal Force Demonstrator

**Product Code :**  
LBNY-0004-10100015



**Description :**

The function of the instrument is electric rotation. The speed is adjustable. The rotating speed and centripetal force can be displayed respectively,  
With memory and zeroing function,  
The weight of the weight can be replaced and the spacing distance is adjustable.  
The product is composed of a base, hand ring track, aluminum frame, steel ball, weights etc.

**Technical Specification :**

The base  $\varnothing$ 235.5mm, height 90.5mm  
Size of suspension ring: 34x12x81mm;  
Aluminum track length 500mm;



**Laboratory Instrument India**