



Email : sales@laboratoryinstrumentindia.com

Product Name :
Hand-held Direct Vision Spectroscope

Product Code :
LBNY-0004-10100016



Description :

The product is composed of an outer cylinder, a telescopic cylinder, prism, slit eyepiece system, etc. Outer cylinder, a telescopic cylinder is made of engineering plastic, the outer cylinder dimensions for $\varnothing 20 \times 60$ mm, wall thickness 1.5mm; The eyepiece system consists of a lens holder, an eyeglass lens seat, size $\varnothing 24 \times 13$ mm; slit width 1.4mm, width 7mm.

Technical Specification :

Product size $\varnothing 20 \times 76$ mm.



Laboratory Instrument India