



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Micro Current Amplifier

**Product Code :**  
LBNY-0004-10100017



**Description :**

The product is composed of the shell, the input end, the output end, the amplification regulation, the power switch and the indicator lamp.  
The product shell is made of ABS engineering plastics, and the input end and the output end adopt plug-in dual-purpose terminal post.

**Technical Specification :**

Instrument use voltage: DC6V; magnification: not less than 600 times; input current: 0.5 ~ 10uA; instrument structure is compact, the experiment is convenient.  
Overall dimensions: 130 \* 60 \* 90mm;



**Laboratory Instrument India**