



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Cathode ray tube (straight tube)

**Product Code :**  
LBNY-0004-10100018



**Description :**

Product is composed of a glass shell, a baffle plate, a bracket, a cathode, an anode and the base; Bezel diameter is not less than 80 mm and the length is not less than 235mm. Tube vacuum range of 0.13 PA ( $1 \times 10^{-3}$  torr) ~ 0.07 PA ( $5 \times 10^{-4}$  torr), the cathode ray tube electrode is made of pure aluminum with no less than 99.7%. Baffle plate is made by a thickness of not less than 0.8 mm, the length of not less than 45mm cross type aluminum plate, flexible and flip to upright on the bracket.

**Technical Specification :**

Cathode ray tube (straight tube)



**Laboratory Instrument India**