



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Resistance Law Experiment Device

**Product Code :**  
LBNY-0004-10100023



**Description :**

The product is composed of a bottom plate, three different kinds of metal wire (copper wire, iron wire, nickel chrome wire), terminal, connecting piece etc.

**Technical Specification :**

Technical parameters of metal wire: a copper wire, diameter:  $0.3\text{mm} \pm 0.04\text{mm}$ , effective length:  $400\text{mm} \pm 2\text{mm}$ ; a iron wire, diameter:  $0.3\text{mm} \pm 0.04\text{mm}$ , effective length:  $400\text{mm} \pm 2\text{mm}$ ; two nickel chrome wire, diameter:  $0.3\text{mm} \pm 0.04\text{mm}$ , effective length:  $400\text{mm} \pm 2\text{mm}$ .  
Base:  $460 \times 150 \times 16\text{mm}$



**Laboratory Instrument India**