



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Digital Timer Air Cushion Guide

**Product Code :**  
LBNY-0004-10100024



**Description :**

The product is a set of mechanical demonstration experiment, the time is measured by digital timer, air cushion guide is used for small gas source, and the products can be used in many kinds of mechanical experiments such as average speed, acceleration and elastic collision.

**Technical Specification :**

Digital Timer Air Cushion Guide



**Laboratory Instrument India**