



Email : sales@laboratoryinstrumentindia.com

Product Name :
Single Turn Coil Motor Principle Demonstration

Product Code :
LBNY-0004-10100044



Description :

The instrument can be used in the experiment of DC motor and AC / DC power generation experiment. In the motor experiment, change the current flow in the coil or the direction of the magnetic field to change the direction of rotation of the coil.

Technical Specification :

At the same time, the red and blue LED lights on the sides of the coil, respectively indicate flowing in and out of the current direction of the coil, the coil is made of 2mm copper tube, coil frame size: 74x100 (±2mm).2, with the J0401 type demonstration meter, 100mA full scale, can do the generator experiment: (rotating coil to do DC power generation experiment, rotating magnetic field to do AC power generation experiment). 3, size: 180mmx235mm.



Laboratory Instrument India