



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Motion Simulation Demonstration Device

**Product Code :**  
LBNY-0004-10100048



**Description :**

The product working voltage AC220V ,composed of a LED simulation screen, an acceleration knob and a deflection knob. Adjust the deflection knob, and the particles do the horizontal projectile motion. By adjusting the acceleration knob, the particles do uniform acceleration of linear motion.

**Technical Specification :**

Motion Simulation Demonstration Device



**Laboratory Instrument India**