



Email : sales@laboratoryinstrumentindia.com

Product Name :
Pointer Electroscop

Product Code :
LBNY-0004-10100058



Description :

This product is composed of two identical pointer electroscopes.
The pointer is made of a thin sheet of metal. The length is not less than 100mm, and the lower part is arrowhead shaped with red lacquer.
The surface diameter is 170mm, and the product support frame is 6mm in diameter, and the iron surface is chrome plated.
The pointer electroscopes are composed of a casing, ball or disc, conductive rod, insulator, pointer, pointer frame, grounding column etc.

Technical Specification :

Pointer Electroscop



Laboratory Instrument India