



Email : sales@laboratoryinstrumentindia.com

Product Name :
Heatconduction

Product Code :
LBNY-0004-10100060



Description :

The product is made up of test stick, heating block and hand held supporting rod. Test bars including: brass rods, rods copper, iron, aluminum bar, stainless steel rods, rod test size is $\varnothing 5 \times 80 \text{mm}$ handheld rod length of 240mm, using a metal rod production, heating block is made of copper material; handheld rod attached to the handle, and the handle adopts insulation materials production, grip length 80mm.

Technical Specification :

Heatconduction



Laboratory Instrument India