



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Resonance Demonstrator

**Product Code :**  
LBNY-0004-10100061



**Description :**

Vibration frequency can be adjusted, vibration frequency can be locked, digital display.  
Motor drive, three spring vibration.  
Using voltage: AC220V + 10%, 50HZ.

**Technical Specification :**

Instrument size: 32CMx15.5CMx56CM.



**Laboratory Instrument India**