



Email : sales@laboratoryinstrumentindia.com

Product Name :
Longitudinal Wave Demonstrator

Product Code :
LBNY-0004-10100063



Description :

Products as a supporting rod suspension spring type, mainly composed of a rack, a spring steel sheet, coil spring, connecting rod, interlining, the source of vibration. The frequencies of spring steel sheet can be in 2.5 ~ 3.3Hz within the scope of adjustment. Frame consists of a left, a right bracket and rubber feet, is made of engineering plastic, the bottom is about 240mm, 25mm wide, frame with a total height of about 380 mm; spring steel sheet thickness of 0.8mm, 300 mm long and 3 mm wide, by changing the oscillator position.

Technical Specification :

Longitudinal Wave Demonstrator



Laboratory Instrument India