



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Micro Deformation Demonstrator

**Product Code :**  
LBNY-0004-10100070



**Description :**

The product is composed of two sets of flat mirrors, a bracket, a laser pointer and a bracket. The output power of the laser pen is less than 1MW, laser pen bracket height is about 35mm, the angle of laser pen on the bracket is adjustable. Plane mirror size is 90 \* 60MM, plane mirror on the bracket height can be adjusted, adjusting range is not less than 40mm.

**Technical Specification :**

Micro Deformation Demonstrator



**Laboratory Instrument India**