



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Pressure and Intensity of Pressure Demonstrator

**Product Code :**  
LBNY-0004-10100071



**Description :**

The instrument is composed of a transparent cup, sponge, pressure block, weight tray, a bearing rod, rulers, etc., instrument appearance size: 79mm×160mm, pressure block for block 3 pcs of 60mm, 40mm, 20mm different diameter which can be replaced, sponge thickness is not less than 38 mm, weight tray diameter 40mm, bearing the pressure rod is made of metal material, Diameter 10mm, effective length is not less than 80mm, scale range 0-6cm.

**Technical Specification :**

Pressure and Intensity of Pressure Demonstrator



**Laboratory Instrument India**