



Email : sales@laboratoryinstrumentindia.com

Product Name :
Demonstration of Micro Current Resistance Meter

Product Code :
LBNY-0004-10100077



Description :

Micro DC detection: DC(G)-50uA-0-500uA 0-100uA. DC voltage:DC(V)0-1V-2.5V-5V-10V-25V-50V-100V-150V-250V resistance measurement:R \times 1??R \times 10??R \times 100??R \times 1K? meter precision \pm 2.5%.

Technical Specification :

Demonstration of Micro Current Resistance Meter



Laboratory Instrument India