



Email : sales@laboratoryinstrumentindia.com

Product Name :
Insulation Resistance Meter

Product Code :
LBNY-0004-10100079



Description :

Swing handle rated speed: 120r/min. The accuracy level is 10 (equivalent to 1 levels of arc length). Product material: Metal. The product is a crank insulation resistance meter. 500V, 25 - 3 type. Rated voltage: 500V + 10%, measured resistance range: 0 ~ 500M.

Technical Specification :

Insulation Resistance Meter



Laboratory Instrument India