



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Insulation Resistance Meter

**Product Code :**  
LBNY-0004-10100079



**Description :**

Swing handle rated speed: 120r/min. The accuracy level is 10 (equivalent to 1 levels of arc length). Product material: Metal. The product is a crank insulation resistance meter. 500V, 25 - 3 type. Rated voltage: 500V + 10%, measured resistance range: 0 ~ 500M.

**Technical Specification :**

Insulation Resistance Meter



**Laboratory Instrument India**