




Email : sales@laboratoryinstrumentindia.com

Product Name :
Demonstration of The Thermometer

Product Code :
LBNY-0004-10100083



www.laboratoryinstrumentindia.com

Description :

The product is composed of a glass thermometer, a scale brand and a protective sleeve. Temperature sensitive liquid: red, packed in a warm bulb and a capillary tube. Scale plate: made of thick 0.5mm aluminum plate, scale plate length 500mm, the scale of the line has two kinds of temperature indication, Celsius temperature (symbol C) and thermodynamic temperature (symbol K). Warm bubble: About 135mm, located in the lower end of the thermometer. Capillary: installed in the scale board, the capillary cross-section of the triangle, when the positive observation from the front, there is amplification. Scale range: 0 ~ 100 degrees Celsius, thermodynamic scale of 273 ~ 373K. Each 10 degree Celsius line length is 4mm. Each 50 degree Celsius for the short term, each 100 degree Celsius for the long term. Showing the value of tolerance ± 1 degree Celsius.

Technical Specification :

Demonstration of The Thermometer



Laboratory Instrument India