



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Electronic Stopwatch

**Product Code :**  
LBNY-0004-10100085



**Description :**

Stopwatch with simple time, split time, and must be with pause button daily timing alarm and automatic snooze function, shows the time, 12 and 24 hour format, calendar, week, shockproof structure function. Using electronic chip, the battery voltage is 1.5V. Specifications for the 110mm x 90mm x 35mm, the data can be accurate to 0.01s.

**Technical Specification :**

Electronic Stopwatch



**Laboratory Instrument India**