



Email : sales@laboratoryinstrumentindia.com

Product Name :
Electronic Stopwatch

Product Code :
LBNY-0004-10100085



Description :

Stopwatch with simple time, split time, and must be with pause button daily timing alarm and automatic snooze function, shows the time, 12 and 24 hour format, calendar, week, shockproof structure function. Using electronic chip, the battery voltage is 1.5V. Specifications for the 110mm x 90mm x 35mm, the data can be accurate to 0.01s.

Technical Specification :

Electronic Stopwatch



Laboratory Instrument India