



Email : sales@laboratoryinstrumentindia.com

Product Name :
Multifunctional Experimental Support

Product Code :
LBNY-0004-1010092



Description :

Products of the general support for teaching, a can to be assembled to form vertical, horizontal and inclined, hanging, clamping, fastening, placement of a variety of experimental support. Various parts of the product are equipped with center sleeve hole M6 * 26mm locking screw. Products are mainly consists of the following components constitute: in the middle of two V-shaped base is connected with a metal rod, 2 pcs vertical rod, 4 pieces bosshead, a universal clamp, the big and the small ring each 1, a flask holder, a short rod, 4 pieces hook roll, floor, a titration clip, a circular tray, test tube rack plate, funnel frame plate, insulation rod and the edge of the table clamp.

Technical Specification :

Multifunctional Experimental Support



Laboratory Instrument India