



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Staff Gauges Plates and Accessories

**Product Code :**  
LBNY-0002-840004



**Description :**

Staff Gauges Plates and Accessories

**Technical Specification :**

Staff Gauges Plates and Accessories

Length: 1.0 m

Wide: 80 mm,  $\pm 10$  mm

Markings: Centimetre markings with labelled decimetres, supply with separate meter number plates

Industrial level precise measurement, reliable results

Stainless steel construction, durable in use and never deform or damage

Mainly used for thickness measurement of leather, paper, film, metal sheet. etc similar thinner materials

Material: 2 mm steel plate with porcelain coated to resist rust and discoloration.

Layout Clear dark markings on a white background with 6 holes of 6 inches wide

Figure Plate (Gauge tags): Figure plates with white background with 2 holes indicating the elevation of the staff gauge plate is installed every meter for every staff gauge. The number of Figure plates are determined by the number of staff gauges to be constructed and ranges of gauges.



**Laboratory Instrument India**

