



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Smart Phone Hardware Component

**Product Code :**  
LBNY-0002-850001



**Description :**

Smart Phone Hardware Component, Smart Phone Hardware Component Equipments, Smart Phone Hardware Component Tools, Smart Phone Hardware Component Tool Kits, Smart Phone Hardware Component Manufacturers, Smart Phone Hardware Component Suppliers from India, China, Kenya

Lab Tender Suppliers Africa and Lab Equipment Africa Kenya, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Benin, Burundi, Tunisia, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Mozambique, Madagascar, South Africa, Ivory Coast, Cameroon, Niger, Burkina Faso, Mayotte, Sao Tome, Seychelles, South Sudan, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Mauritania, Eritrea, Gambia, Botswana, Namibia, Gabon, Lesotho, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia

**Technical Specification :**

Smart Phone (Hardware Component)  
Smart-Phone: Any standard electronic consumer smart-phone (not proprietary)  
Operating system: Android OS 8.0 or iOS 9+  
Network: GSM 850 / 900 / 1800 / 1900 3G bands HSDPA 850 / 900 / 1900 / 2100 4G  
Internet connectivity: 512 kb/s  
RAM Memory: 2GB  
Data Storage: MicroSDXC 16 GB  
Components: Wi-Fi 802.11 a/b/g/n/ac  
Camera: 13 MP  
GPS sensor A-GPS, GLONASS

---

Settings: Camera, high accuracy GPS location and permissions granted

Battery > 3000 mAh Li-Ions

Display Glass front LCD 4.7" (11.94 cm) dimensions

Resolution: 720 x 1280 pixels

SIM Card: From Available mobile GSM provider

Power supply: Battery charger, cables, and accessories

Protection: Rubber case and screen protector



**Laboratory Instrument India**