



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Object ups and downs condition demonstrator

**Product Code :**  
LBNY-0004-101000112



**Description :**

Liquid containing barrel is a transparent plastic cylinder, the internal diameter of the shell is not less than 95mm, depth is not less than 285mm. The product is composed of a liquid container, a floating body and accessories. Floating body is composed of a diving boat model, buoy, a counterweight body. The floating body can be floating, suspended, or sinking state in a liquid. The submarine size about 100 x 30 x 35mm; float size: 40mm x 100mm. The counterweight body is composed of 4 supporting and weights 20 grams ; accessories include skin tube clamp, density meter, sucking ball etc.

**Technical Specification :**

Object ups and downs condition demonstrator



**Laboratory Instrument India**