




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Experimental box of middle school electricity

**Product Code :**  
LBNY-0004-101000116



**Description :**

It is used to compose series circuit and parallel circuit, measure current intensity with ammeter, measure voltage with voltmeter, change current density with slide rheostat, measure resistance and power of bulb with voltmeter and ammeter and self-combination of experiment.

**Technical Specification :**

Dimension (L\*W\*H CM): 23.2\*15.4\*5.4



**Laboratory Instrument India**