



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Field Computer Hardware Component

**Product Code :**  
LBNY-0002-850002



**Description :**

Field Computer Hardware Component

**Technical Specification :**

Field Computer (Hardware Component)  
Dimensions (mm): 272.00 x 345.00 x 32.80  
Weight (kg): 1.5-2.2  
Operating system: Windows 10 Professional  
Battery Life (up to hours): 20  
Display  
Size: 14.00 - 17.00 inch  
Processor  
11th Generation Intel® Core™ i7-1165G7@ 2.5GHz 2.5 GHz  
Cache: 8MB or better  
Memory  
RAM: 16GB or above  
Storage  
SSD: 1TB  
Connectivity  
Bluetooth version: 5.0 or above  
Inputs  
Web Camera: Yes  
Pointer Device: Touchpad

---

Backlit Keyboard: Yes

Touchpad: Yes

Internal Mic: Yes

Speakers: 2 Speakers

Ports and slots

Number of USB Ports: 4

USB Ports: 1 x USB 2.0, 3 x USB 3.1 Gen 1 (Type A)

HDMI Port: Standard

Multi Card Slot : SDXC Card Reader

Headphone and Mic Combo Jack: Yes

RJ45 (LAN): Yes



**Laboratory Instrument India**