



Email : sales@laboratoryinstrumentindia.com

Product Name :
Plane Lift Demonstrator

Product Code :
LBNY-0004-101000117



Description :

It is used to demonstrate and verify the principle of airplane wing producing lift when the airplane flies in the atmosphere. It consists of fan, airplane model, transparent cover, slide bar, slide bar, clamp plate and board base.

Technical Specification :

Dimension (L*W*H CM): 52*15*27

Inner packs dimension (L*W*H CM): 53.5*18*28.5



Laboratory Instrument India