



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Plane Lift Demonstrator

**Product Code :**  
LBNY-0004-101000117



**Description :**

It is used to demonstrate and verify the principle of airplane wing producing lift when the airplane flies in the atmosphere. It consists of fan, airplane model, transparent cover, slide bar, slide bar, clamp plate and board base.

**Technical Specification :**

Dimension (L\*W\*H CM): 52\*15\*27  
Inner packs dimension (L\*W\*H CM): 53.5\*18\*28.5



**Laboratory Instrument India**