



Email : sales@laboratoryinstrumentindia.com

Product Name :
Fluorescent Lamp Demonstrator

Product Code :
LBNY-0004-101000134



Description :

It is used for student to understand application of self-induction and effect of fluorescent lamp components, so as to grasp installation method of fluorescent lamp. It consists of fluorescent tube, ballast, starter and switch.

Technical Specification :

Dimension (L*W*H CM): 40.8*4*21

Inner packs dimension (L*W*H CM): 42.5*7.9*21.5



Laboratory Instrument India