



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Fluorescent Lamp Demonstrator

**Product Code :**  
LBNY-0004-101000134



**Description :**

It is used for student to understand application of self-induction and effect of fluorescent lamp components, so as to grasp installation method of fluorescent lamp. It consists of fluorescent tube, ballast, starter and switch.

**Technical Specification :**

Dimension (L\*W\*H CM): 40.8\*4\*21

Inner packs dimension (L\*W\*H CM): 42.5\*7.9\*21.5



**Laboratory Instrument India**