



Email : sales@laboratoryinstrumentindia.com

Product Name :
Satellite Transmitter

Product Code :
LBNY-0002-860001



Description :

Satellite Transmitter, Satellite Transmitter Equipments, Satellite Transmitter Tools, Satellite Transmitter Tool Kits, Satellite Transmitter Manufacturers, Satellite Transmitter Suppliers from India, China, Kenya

Lab Tender Suppliers Africa and Lab Equipment Africa Kenya, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Benin, Burundi, Tunisia, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Mozambique, Madagascar, South Africa, Ivory Coast, Cameroon, Niger, Burkina Faso, Mayotte, Sao Tome, Seychelles, South Sudan, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Mauritania, Eritrea, Gambia, Botswana, Namibia, Gabon, Lesotho, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia

Technical Specification :

Satellite Transmitter
Baud rate: 1.200
Frequency Range: 401 to 403 MHz
Frequency Stability: 0.25 PPM
Timing accuracy: +/- 0.5 seconds
Modulation Offset: QPSK
Channel bandwidth: 3KHz
Power: 12.5 Vdc nominal
Power Consumption: 30 mA during GPS use / 3 Amps max at 10 W RF output power
RF Power: 1 – 5 watts (adjustable)
GPS receiver: Embedded

Connectors: GPS antenna, RF output, power supply, serial ports
User interface: RS-232 and/or RS-485 serial port
Software: Programming software for parameters configuration and control



Laboratory Instrument India