



Email : sales@laboratoryinstrumentindia.com

Product Name :
Ultrasonic Application Demonstrator

Product Code :
LBNY-0004-101000137



Description :

It performs demonstration experiment by using professional microprocessor control circuit, by ultrasonic transmitting and receiving and relevant control circuit. It is used for the demonstration of control experiments: ultrasonic ranging, alarm, water level control and car back alarm.

Technical Specification :

Dimension (L*W*H CM): 19.4*17.6*6.9
Inner packs dimension (L*W*H CM): 21.1*21.4*9



Laboratory Instrument India