



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Ultrasonic Application Demonstrator

**Product Code :**  
LBNY-0004-101000137



**Description :**

It performs demonstration experiment by using professional microprocessor control circuit, by ultrasonic transmitting and receiving and relevant control circuit. It is used for the demonstration of control experiments: ultrasonic ranging, alarm, water level control and car back alarm.

**Technical Specification :**

Dimension (L\*W\*H CM): 19.4\*17.6\*6.9  
Inner packs dimension (L\*W\*H CM): 21.1\*21.4\*9



**Laboratory Instrument India**