



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Charger of Various Functions

**Product Code :**  
LBNY-0004-101000140



**Description :**

It is used for the maintenance and automatic-charging of lead-acid battery, adjustable internal resistance battery and common lead-acid battery.  
As well as, it is used for the automatic-charging of rechargeable battery in different trade.

**Technical Specification :**

Charger of Various Functions



**Laboratory Instrument India**